

Colour options



B - brilliant blue



P - orange



F - ocean green



M - corn yellow

Technical Parameters

| | PURA Standard | PURA Comfort/Exclusive |
|--|---------------|------------------------|
| Length of loading area | 2100 mm | 2100 mm |
| Width of loading area | 590 mm | 590 mm |
| Height of loading area, including upholstery | | |
| Minimum | 630 mm | 610 mm |
| Maximum | 630 mm | 1100 mm |
| Trendelenburg position | -12° | -12° |
| Armrest adjustment angle | | |
| Vertical | +15°/-30° | +15°/-30° |
| Horizontal | -10°/-170° | -10°/-170° |
| Maximum patient weight | 180 kg | 180 kg |


Types of PURA dialysis chair

DIA-AX – PURA Standard

DIA-BX – PURA Comfort

DIA-CX – PURA Exclusive

X – letter indicating upholstery colour



DIALYSIS CHAIR



BORCAD cz s.r.o.

CZ – 739 45 Fryčovice 673

Tel.: +420 558 640 631

Fax: +420 558 668 087

e-mail: borcad@borcad.cz

www.borcad.cz

BORCAD®

DIALYSIS CHAIR




Quick cleaning thanks to seamless upholstery.

Pleasant design combined with compact proportions for easy treatment.

Very comfortable loading area makes the patient's stay more pleasant.

Wide ergonomically adjustable armrest increases the patient's comfort.



BORCAD®

Dialysis chair PURA

The attractive design of PURA creates a feeling of wellbeing for the patient.

This unique product design sets high standards and increases the prestige of a dialysis centre.

Designed to take up the minimum amount of floor space to maximize treatment areas.

Unique design combined with a high degree of comfort.




Wide ergonomically adjustable armrest increases the patient's comfort and independence. Armrest could be controlled by client which makes the work of hospital staff easier.



Removable armrests for transporting the chair through narrow areas.



Unique foot support contributes to the patients comfort and allows the patient to change position with the minimal amount of effort.



Compact seamless areas speed up cleaning and enable disinfection with no compromise.












Unambiguous and easy control without time consuming training.

The hand control keyboard can be locked using magnetic key.



Versions of the PURA chair

| | PURA-Standard (DIA-AX) | PURA-Comfort (DIA-BX) | PURA-Exclusive (DIA-CX) |
|----------------------------|--|--|--|
| Height adjustment |  Fixed height |  Adjustable height |  Adjustable height |
| Leg rest leaning | Leg rest fixed connected with the sitting part | Leg rest leaned independently of the sitting part | Leg rest leaned independently of the sitting part |
| Castor versions |  Average of 75 mm |  Average of 100 mm |  Average of 100 mm |
| Foot support controlled by |  staff |  staff as well as patient |  staff as well as patient |
| Chair fixation | Each castor is being braked separately | Each castor is being braked separately | All castors are braked by the central brake |

Accessories



DIA-ZP-009.X Adjusting foot support



DIA-ZP-004.X Head rest



DIA-ZP-001.X Head rest, semisecular



DIA-ZP-002 Handrail



DIA-ZP-003 Infusion stand holder



SLA-ZP-003 Telescopic infusion stand



DIA-ZP-006 Lamp supplied along with handrail



DIA-ZP-010 PVC cover for foot support